```
In [2]: %matplotlib inline
import matplotlib.pyplot as plt

In [19]: ##starting up the SURFv5...

##
## The board has just been programmed
##
## The calibration file is empty:
## the cal file calibrations/surf_calibrations.json is an
## empty file with just 2 curly brackets:
## {}
## import the useful functions:
import surf_board
import calibrations.surf_calibrations as sc
reload(sc)

Out[19]: <module 'calibrations.surf_calibrations' from 'calibrations/surf_calibrations.py'>

In [10]: ## verify that the surf_calibrations.json file is, in fact, empty:
sc.read_cal() #when calling this function without any arguments,
#it returns the name of the boards in the cal file
```

```
In [11]: ## run the surf board module for the board 'CANOES':
         surf_board.do('CANOES')
         Identification Register: 53354137 (S5A7)
         Version Register: 0.2.12 compiled 22/27
         Device DNA: 70b0cc13d9705c
         path: /sys/class/uio/uio0
         Using default VadjN
         Using default VadjP of 2700.
         setting default values
         Found sync edge: 2469
         Found WR_STRB edge on LAB0: 3744
         Found WR_STRB edge on LAB1: 3770
         Found WR_STRB edge on LAB2: 3756
         Found WR_STRB edge on LAB3: 3767
         Found WR_STRB edge on LAB4: 3738
         Found WR_STRB edge on LAB5: 3800
         Found WR_STRB edge on LAB6: 3820
         Found WR_STRB edge on LAB7: 3822
         Found WR_STRB edge on LAB8: 3720
         Found WR_STRB edge on LAB9: 3696
         Found WR_STRB edge on LAB10: 3810
         Found WR_STRB edge on LAB11: 3906
         adding board to calibration json file..
         scanning for Vadjp..
         LAB4D# 0 Trial: vadjp 2720 width 602.000000 target 543.000000
         LAB4D# 0 Trial: vadjp 2710 width 565.000000 target 543.000000
         LAB4D# 0 Trial: vadjp 2700 width 522.000000 target 543.000000
         LAB4D# 0 Trial: vadjp 2705 width 543.000000 target 543.000000
         LAB4D# 1 Trial: vadjp 2680 width 511.000000 target 542.000000
         LAB4D# 1 Trial: vadjp 2690 width 550.000000 target 542.000000
         LAB4D# 1 Trial: vadjp 2685 width 532.000000 target 542.000000
         LAB4D# 1 Trial: vadjp 2687 width 540.000000 target 542.000000
         LAB4D# 2 Trial: vadjp 2720 width 561.000000 target 531.000000
         LAB4D# 2 Trial: vadjp 2710 width 516.000000 target 531.000000
         LAB4D# 2 Trial: vadjp 2715 width 538.000000 target 531.000000
         LAB4D# 2 Trial: vadjp 2713 width 531.000000 target 531.000000
         LAB4D# 3 Trial: vadjp 2720 width 604.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2710 width 566.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2700 width 526.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2705 width 545.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2703 width 541.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2701 width 530.000000 target 535.000000
         LAB4D# 3 Trial: vadjp 2702 width 534.000000 target 535.000000
         LAB4D# 4 Trial: vadjp 2680 width 483.000000 target 556.000000
         LAB4D# 4 Trial: vadjp 2690 width 531.000000 target 556.000000
         LAB4D# 4 Trial: vadjp 2700 width 567.000000 target 556.000000
         LAB4D# 4 Trial: vadjp 2695 width 548.000000 target 556.000000
         LAB4D# 4 Trial: vadjp 2697 width 556.000000 target 556.000000
         LAB4D# 5 Trial: vadjp 2720 width 507.000000 target 507.000000
         LAB4D# 6 Trial: vadjp 2720 width 518.000000 target 494.000000
         LAB4D# 6 Trial: vadjp 2710 width 479.000000 target 494.000000
         LAB4D# 6 Trial: vadjp 2715 width 497.000000 target 494.000000
         LAB4D# 6 Trial: vadjp 2713 width 488.000000 target 494.000000
         LAB4D# 6 Trial: vadjp 2714 width 493.000000 target 494.000000
         LAB4D# 7 Trial: vadjp 2680 width 454.000000 target 509.000000
         LAB4D# 7 Trial: vadjp 2690 width 496.000000 target 509.000000
         LAB4D# 7 Trial: vadjp 2700 width 519.000000 target 509.000000
         LAB4D# 7 Trial: vadjp 2695 width 506.000000 target 509.000000
         LAB4D# 7 Trial: vadjp 2697 width 512.000000 target 509.000000
         LAB4D# 7 Trial: vadjp 2696 width 507.000000 target 509.000000
LAB4D# 8 Trial: vadjp 2720 width 624.000000 target 584.000000
         LAB4D# 8 Trial: vadjp 2710 width 580.000000 target 584.000000
         LAB4D# 8 Trial: vadjp 2715 width 603.000000 target 584.000000
         LAB4D# 8 Trial: vadjp 2713 width 591.000000 target 584.000000
         LAB4D# 8 Trial: vadjp 2711 width 585.000000 target 584.000000
         LAB4D# 9 Trial: vadjp 2720 width 669.000000 target 614.000000
         LAB4D# 9 Trial: vadjp 2710 width 619.000000 target 614.000000
         LAB4D# 9 Trial: vadjp 2700 width 579.000000 target 614.000000
         LAB4D# 9 Trial: vadjp 2705 width 599.000000 target 614.000000
         LAB4D# 9 Trial: vadjp 2710 width 621.000000 target 614.000000
         LAB4D# 9 Trial: vadjp 2708 width 613.000000 target 614.000000
         LAB4D# 10 Trial: vadjp 2720 width 576.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2710 width 545.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2700 width 506.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2705 width 527.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2703 width 517.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2704 width 518.000000 target 521.000000
         LAB4D# 10 Trial: vadjp 2705 width 523.000000 target 521.000000
         LAB4D# 11 Trial: vadjp 2720 width 509.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2710 width 460.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2700 width 417.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2705 width 440.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2710 width 461.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2708 width 452.000000 target 447.000000
         LAB4D# 11 Trial: vadjp 2706 width 445.000000 target 447.000000
         Saving VadjP for board CANOES
         no entry matching input key
          setting default values
         Clock Status: LAB4 Clock is enabled (CLKSEL[1] = 1)
                      : LAB4 Driving Clock is FPGA Clock (CLKSEL[0] = 0)
                      : Local Clock is enabled (CLKSEL[2] = 0)
                      : FPGA System Clock PLL is running (PLLCTRL[0] = 0/PLLCTRL[2] = 0)
                      : FPGA System Clock is Local Clock (PLLCTRL[1] = 0)
          Int Status : 00000000
                     : Internal value 000, Key value 000
          I FD
                     : 00000000
          Full LED
                      : 00000000
          Int Mask
         Reading from MCP4728...
         DAC channel A (RFP VPED 0): register is set to 0x9c4, EEPROM is set to 0x9c4
         DAC channel B (RFP_VPED_1): register is set to 0x578, EEPROM is set to 0x578 DAC channel C (RFP_VPED_2): register is set to 0x578, EEPROM is set to 0x578
                                    : register is set to 0x76c, EEPROM is set to 0x76c
         DAC channel D (VPED)
         LAB4 runmode: not enabled
         LAB4 testpat: not enabled
Out[11]: <SURF at /sys/class/uio/uio0>
In [13]: | ## recheck the surf_calibrations.json file
         sc.read_cal()
```

board in cal file: CANOES

```
In [15]: ## load up the surf_data class:
           import surf_data
           reload(surf_data)
           dev=surf_data.SurfData()
In [20]: ## check the calibration data available for CANOES:
          sc.read cal(dev.dev.dna())
          vadjp
DNA
In [21]: | ## the board DNA, and vadjP are in the cal file
           ## let's take some pedestal data:
           dev.pedestalRun()
          Saving pedestals for board CANOES...
Out[21]: array([[1638, 1809, 1631, ..., 1719, 1654, 1660],
                  [1645, 1828, 1633, ..., 1748, 1665, 1711],
[1632, 1827, 1607, ..., 1721, 1664, 1699],
                   [1640, 1874, 1612, ..., 1733, 1693, 1497],
                  [1616, 1860, 1637, ..., 1740, 1664, 1467],
[1633, 1837, 1623, ..., 1698, 1669, 1454]])
In [22]: ##OK, lets check the cal file again for CANOES:
           sc.read_cal(dev.dev.dna())
          pedestals
          vadjp
DNA
In [23]: ## good, the pedestal data was added as expected
          ## now let's check to see if the data is sensible. Take 5 events:
data = dev.log(5, save=False)
In [30]: for i in range(11):
          plt.plot(data[1][i], 'o--', ms=2, label='CH{:d}'.format(i))
plt.xlabel('event sample')
           plt.ylabel('amplitude [lab4d counts]')
           plt.ylim([-30, 50])
           plt.legend(ncol=3)
Out[30]: <matplotlib.legend.Legend at 0xbf30aec>
                50
                                     - - CH0
                                                   - → CH4
                                                                 ← → CH8
                40
                                     - → CH1

← → CH5

                                                                 - → CH9
                                                                 - - CH10
                                     - → CH2
                                                   - - CH6
           amplitude [lab4d counts]
                30
                                                   - - CH7
                                     - → CH3
                20
```

In []:

10

-10

-20

-30 L

200

400

600

event sample

1000

1200

800